

## Image Guided Radiosurgery for Trigeminal Neuralgia

Ning Wen<sup>1</sup>, Yimei Huang<sup>2</sup>, Karen Snyder<sup>1</sup>, Haisen Li<sup>3</sup>, Yujiao Qin<sup>4</sup>, Essa Mayyas<sup>2</sup>, Indrin J. Chetty<sup>1</sup>, M. Salim Siddiqui

1. Radiation Oncology, Henry Ford Health System, Detroit, USA 2. Medical Physics, Henry Ford Health System 3. Medical Physics, Henry Ford Health 4. Medical Physics, Advocate Health Care, Chicago, Select Country 5.

✉ **Corresponding author:** Ning Wen, nwen1@hfhs.org

**Categories:** Medical Physics

**Keywords:** stereotactic radiosurgery, trigeminal neuralgia, linac, cone beam ct, edge

**How to cite this abstract**

Wen N, Huang Y, Snyder K, et al. (June 16, 2016) Image Guided Radiosurgery for Trigeminal Neuralgia. Cureus 8(6): a119

## Abstract

**Objectives:** To share our initial clinical experiences for the treatment of trigeminal neuralgia (TN) on a novel image-guided stereotactic linear accelerator.

**Methods:** The latest LINAC-based SRS platform (the Edge™, Varian Medical Systems, Palo Alto, CA) has features including high intensity mode (6 MV flattening filter free mode with a maximum dose rate of 1400 MU/minute), a head frame system, a 6 DoF robotic couch, and conical cones. The diameters of the conical cones range from 4 to 17.5 mm. The system also has an advanced imaging package for treatment localization, including Cone Beam Computed Tomography (CBCT) for volumetric assessment, stereoscopic x-ray imaging and an optical surface monitoring system (OSMS) for monitoring patients via 3D surface imaging. The end-to-end hidden target test was performed using a Rando head phantom that was attached to the treatment table via a head frame and localized with CBCT. The planar dose distribution from arc delivery with a 4mm cone was verified using the Lucy phantom (Standard Imaging) with Gafchromic EBT3 films. The patient was simulated in a Philips Brilliance big bore CT scanner (1 mm slice thickness; Matrix 512x512, 120kVp, 501mAs). A high resolution MR simulation (BFFE, 0.5 mm slice thickness) was also acquired for contouring purpose. The treatment plan was done in the Eclipse cone planning system. 14 partial arcs with 4 mm cone were included in the treatment plan for a dose of 85-90 Gy at isocenter. In the treatment room, the patient was first aligned with OSMS. CBCT (full fan, full rotation, 1 mm slice thickness) was used for patient localization and 6 DoF corrections were applied subsequently. A second CBCT was acquired afterwards for verification. Intra-fraction CBCT was acquired mid-treatment at couch zero to evaluate intra-fractional motion. The OSMS was used to monitor patient motion throughout the treatment with 1mm/1° tolerance.

**Results:** The overall system accuracy using 15 mm cone was 0.34±0.11 mm based on End-to-End test. The 4 mm cone was concentric with the 15 mm cone within 0.15 mm. The output of a single 160° arc in 4 mm cone measured with EBT3 film was consistent with that from the initial commissioning measurement. Two patients with TN have been treated and both experienced favorable outcome. Pre-treatment verification CBCT matched well with reference planning CT within 0.2 mm/1°. The patient movement was within 0.5 mm/1° based on the intra-fraction CBCT and under 1mm/1° from real time surface imaging.

**Conclusions:** CBCT imaging coupled with 6 DoF robotic couch has been shown to be accurate

Open Access

Abstract

Published 06/16/2016

Copyright

© Copyright 2016

Wen et al. This is an open access article distributed under the terms of the Creative Commons Attribution License CC-BY 3.0., which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

Distributed under

Creative Commons CC-BY 3.0

and reproducible through end-to-end phantom measurements. Frame-based patient immobilization minimize intra-fraction motion as demonstrated by the intra-fraction CBCT from the patient treatments. Therefore, the Edge platform combines the gold-standard of patient immobilization and the most advanced online image-guidance. Our initial experience has shown that Edge radiosurgical system is effective for treating TN patients.